

SimpleScore User Manual

1. Installation

Open **MATLAB** (recommended version **R2021a or later**).


Open the **APP** tab in the MATLAB menu bar and select **Install App**.

Choose the provided file **simplescore.mlappinstall** to install the application.

After installation is complete, the app



will appear under **My Apps** in the APP menu.




2. Usage

2.1 Load Data

Open the software and click **[Open File]** to select the data file to be analyzed.

During file loading, the button will appear **dark gray**, indicating that the operation is in progress.



名称	修改日期	类型	大小
20170308-KO-D.btn	2017/3/23 21:04	BTN 文件	9,804,414...

名称: 20170308-KO-D.btn
 类型: (*.btn)
 打开(O) 取消

If no file is selected or if the file format is incorrect, the message “**No File**” will be displayed




2.2 Channel Selection


Select the **EEG** and **EMG** channels from the drop-down menus.

Click **[Read Data]** to load the selected channels.

During data reading, the button will appear **dark gray**.

Once reading is complete, the message “**read over**” will be displayed.






2.3 Feature Calculation

Select the **time window size**, then click [**Calculate Features**] to segment the data and compute feature values.

During computation, the button will appear **dark gray**.

After completion, the calculated features will be displayed and plotted automatically.




2.4 Sleep Stage Scoring

Select [WAKE] / [SWS] / [REM] for manual sleep stage scoring.


- **Left mouse click:** select the boundary of a stage
- **Right mouse click:** assign the sleep stage label

If the manual scoring result is unsatisfactory, click [ClearScore] to clear the current classification and re-score.



2.5 Sleep Stage Visualization

Click [PlotAll] to plot sleep staging results for the **entire recording duration**.



2.6 Sleep Stage Fine Adjustment


If certain staging results are unsatisfactory, select **fine adjustment** mode.

Click [PlotSingle] to examine:

- EEG waveform
- EEG power spectral density (EEG-PSD)
- EMG waveform
- EMG RMS (red line indicates the current epoch)

Use:

- [Previous] / [Next] to navigate between epochs
- [W] / [S] / [R] to change the sleep state




2.7 Analysis and Visualization of Results

Select the **bin range**, then click **[Figure All]** to generate:

- Spectrograms
- Sleep score plots
- EEG plots
- EMG plots


for the selected time range.



Click **Sleep Result Statistics** to obtain summary statistics.

The output includes:


1. A bar plot showing the total duration (in seconds) of **WAKE**, **SWS**, and **REM** within the selected time range
2. A plot showing the number of transitions:
 - WAKE → SWS
 - SWS → REM
 - REM → SWS
 - SWS → WAKE



Bin from To

Time(second) from To

Click **PSD Statistics** to obtain power statistics of different frequency bands under **WAKE**, **SWS**, and **REM** states.



2.8 Saving Sleep Scoring Results

Click [Save] to save the results.


The results are saved in the same folder as the original data, generating two files:

- **score.mat**
Contains the original data, extracted features, and sleep scoring results, facilitating future extensions of scoring functionality.
- **score.xlsx**
Contains the exported sleep scoring results in spreadsheet format.


data > 20170308-KO-D			
搜索"20170308-KO-D"			
名称	修改日期	类型	大小
score.mat	2021/12/7 9:52	MAT 文件	980,743 KB
score.xlsx	2021/12/7 9:51	Microsoft Excel ...	61 KB

2.9 Merging BTN Files



Create an **Excel file** and list, in order, the names of the data files to be merged.



Click **Merge BTN File** and select the previously created Excel file.




Select the corresponding **EEG** and **EMG** channels.




Click **Save Merged** to save the merged file.

⚠ *Note: When the data size is large, saving may take some time.*



A file named **merge.mat** will be saved in the original data folder.

Click **Open MAT File** to open the merged



MAT-format file.